Heat as a form of energy Converting electrical energy into heat

Converting electrical energy into heat energy - Measuring with a voltmeter and an ammeter

## Objects of the experiment

■ The aim of the experiment is to establish the equivalence of electrical and heat energy.

## **Theory**

Water will be heated up and from the temperature difference we will calculate the heat energy. The electrical power will be measured with two multimeters and with the duration we can calculate the electrical energy.

We will use a Dewar vessel to minimize the losses of heat during the experiment. It must additionally be taken into account that the vessel itself possesse a certain heat capacity  $c_D$  (dependent on the liquid height). In order to determine  $c_D$  water with a mass  $m_2$  is poured into the calorimeter and the initial temperature  $\vartheta_2$  is measured after heat compensation.

A quantity of water  $m_1$  with higher temperature  $\vartheta_1$  is then added. The mixture temperature  $\vartheta_m$  is again obtained after heat compensation. The quantity of heat given off is described by

$$\Delta Q_1 = c_1 m_1 (\upsilon_1 - \upsilon_m), \tag{I}$$

while the absorbed heat quantity is

$$\Delta Q_2 = \left(c_D + c_2 m_2\right) \left(\upsilon_m - \upsilon_2\right). \tag{II}$$

Where  $c_1 = c_2 = c_W$  (specific heat capacity of water), we obtain the following :

$$K = \frac{c_D}{c_W} = \frac{m_1(\nu_1 - \nu_m) - m_2(\nu_m - \nu_2)}{(\nu_m - \nu_2)}$$
 (III)

The quotient K can be seen as the mass of the water quantity which has the same heat capacity as the calorimeter.

The electrical energy can be calculated by

$$E = UI \cdot \Delta t \tag{IV}$$

with the voltage U the current I and the time difference.

If the temperature increases from  $\vartheta_1$  to  $\vartheta_2$ , we can set :

$$E = UI \cdot \Delta t = c_w \left( K + m_w \right) \left( \upsilon_2 - \upsilon_1 \right) \tag{V}$$

| Apparatus  1 Electric calorimeter attachement 384 20        |   |
|---|---|
| 1 Dewar vessel  |   |
| Version (a)   |   |
| 1 Thermometer, -10 to 110 C                                 |   |
| 1 Digital thermometer with 1 input                          |   |
| 1 Temperature sensor NiCr-Ni 666 19                         | 3 |
| or - Version (c)  | _ |
| 1 Mobile-CASSY  |   |
| 1 NiCr-Ni temperature sensor 1.5 mm 529 67                  |   |
| 1 Stopclock313 07   |   |
| 1 Beaker, 250 ml squat                                      |   |
| 1 Graduated cylinder, 250 ml 665 75 1 Multimeter LDanalog20 |   |
| 1 Multimeter LDanalog 30531 13                              |   |
| 1 Variable extra low-voltage transformer 521 35             |   |
| 3 Connecting lead 50 cm black501 28                         |   |
| 1 Pair cables 50 cm red/blue 501 45                         |   |

## Measuring example

In a mixing experiment we can measure the water equivalent of the calorimeter. We will get a value of around 9.5 g. Please note that this value is a function of the filling of the calorimeter. Use always the same filling level or make a measurement as a function of the level.

| Mass of the water   | $m_{\text{W}}$ | = 170 g    |
|---------------------|----------------|------------|
| Initial temperature | $\vartheta_1$  | = 20.05 °C |
| Final temperature   | $\vartheta_2$  | = 25.10 °C |
| Time difference     | t              | = 300 s    |
| Current             | I              | = 2.72 A   |
| Voltage             | U              | = 4.7 V    |

## **Evaluation and results**

We can calculte the electrical energy  $E_{\text{elec}} = U \ I \ t = 3835 \ Ws = 3835 \ J$ . With a literature value of  $c_W = 4.2 \ J/gK$  for the specific heat of water we get the heat energy  $E_{\text{heat}} = 3807 \ J$ .




Fig. 1: Experimental setup for measuring with the thermometer 382 34 (Version (a))




Fig. 2: Experimental setup for measuring with the digital thermometer 666 190 (Version (b))




Fig. 3: Experimental setup for measuring with Mobile CASSY 524 009 (Version (c))