

A close-up photograph of various LED lighting components. In the upper left, a circular ring of numerous small, bright white LEDs is visible. To its right, a flexible black cable with a braided outer sheath leads to a cylindrical black light fixture that is emitting a bright white beam of light. Below these, a large, bright, circular white light source is partially visible. The background is a dark, solid blue.

LED ILLUMINATION



MODULAR LED LIGHTING SYSTEM

PHOTONIC develops LED lighting systems for numerous applications with the constant aim of offering customers „up to date“ technologies.

Photonic secures the success of its customers with innovative ideas. The comprehensive modular system for microscopy covers light sources for just about every task and is compatible with all current microscopes. The interaction between intelligent solutions and high product quality is of particular importance to Photonic as this achieves a balanced price-performance ratio.

The LED product line for direct object illumination was developed specifically by Photonic for stereomicroscopy illumination tasks. The modularity of this system enables an upgradable and cost-effective light solution for the following application areas:

- ▶ Automotive industry
- ▶ Industry
- ▶ Forensics
- ▶ Schools and universities
- ▶ Life Sciences
- ▶ Laboratories
- ▶ Measurement and control engineering
- ▶ Quality assurance

LED ILLUMINATION VS. FIBRE OPTICS

One advantage of Photonic LED ringlights compared to fibre optic lighting systems is the possibility of switching segments of the ringlight on/off and to rotate the segments selected in this manner.

In addition to homogeneous illumination, this also makes it possible to achieve directed lighting. Resulting shadow formation also gives rise to excellent contrasting options for surface structures.




The following product properties also offer further advantages compared to fibre optics:

- ▶ Daylight quality (5,000 K)
- ▶ Long service life (up to 25,000 hours)
- ▶ Noise and vibration-free operation
- ▶ Cold light
- ▶ Homogeneous and shadow-free illumination

TYPE AND DIRECTION of lighting are decisive for results. Three approaches are used for the required optimal image contrast:

Incident light (with ringlight [1] or spot [2], offers highest light intensity with maximum flexibility and largest choice of lighting systems) and backlight [3].







LED RINGLIGHTS



LED RINGLIGHT 38/20

This LED ringlight (20 LEDs), with an internal diameter of 38 mm, is suitable for video system and image processing lenses. Using the control unit with segment control provides additional lighting variants.

Working distance: 50–100 mm.



LED RINGLIGHT 66/40

The LED ringlight (40 LEDs) provides exceptionally bright and homogeneous object illumination and is suitable for microscopes with a lens diameter of 66 mm as well as for other microscopes in combination with an adapter. The ringlight is attached with spring and screw, protecting the lens. In combination with the control unit and segment control, various illumination scenarios can be implemented.

Working distance: 55–120 mm.



LED RINGLIGHT 66/80

This LED ringlight has 80 diodes and a correspondingly high light output. Different segments can be selected, lighted and even rotated in combination with the control unit 80. This guarantees optimal detection of contrasts on different surface structures.

Working distance: 35–120 mm or 120–300 mm.

ACCESSORIES

Photonic offers polarizing filter sets as accessories for LED ringlights, comprising an analyzer and a polarizer. The polariser is fastened outside the ringlight and the analyser is integrated in the ringlight. The sets are available for the LED ringlights 66/40 and 66/80.

Polarizing filter set for suppressing unwanted reflections, for ringlights 66/40 and 66/80



Diffuser for light scattering for ringlights 66/40 and 66/80



Adapter for attaching LED ringlights to various microscope lenses (see Accessories on Page 13).







LED HI-POWER SPOTS



LED HI-POWER SPOT

This compact spot is required for all microscopy applications where a small area (15 mm) needs to be illuminated. The very bright light is directly transmitted to the object and the Hi-Power spot is also impressive with a colour temperature of 5,500 K. The light output is comparable to that of a 150 W halogen light source and a fibre-optic cable. One or two Hi-Power spots can be connected to the control unit. If two spots are attached, it is possible to actuate one spot, alternating spots or both spots at the same time.



LED HI-POWER 2-ARM SPOT

Can be attached directly to the microscope stand. Attachment to the stand results in a space-saving application.

- ▶ Cable fed into the flexible arms
- ▶ Neutral white light 5,500 K
- ▶ Homogenous illumination
- ▶ High quality flexible arms
- ▶ Stable mounting on the microscope stand

ACCESSORIES

Diffuser for light scattering,
polarizing filter for suppressing unwanted reflections



Spot holder available in 85, 200 and 300 mm lengths

Base plate for attaching spot holder on which the LED spot lights can be mounted



Adapter for attaching 2-arm spots to columns and various stands (see Accessories on Page 13)





LED SPOTS



LED SPOT 19

The LED spot (19 LEDs) can be directly attached to the microscope with a spot holder or positioned anywhere with a spot holder and base plate. Two LED spots can be connected to and controlled by a single control unit.



LED SPOT 19 2-ARM

This 2-arm LED spot (19 LEDs each) consists of two connected spot holders, each with an LED spot 19 attached to the end. The lighting system can also be directly attached to the microscope stand with an adapter.

ACCESSORIES

Colour filters available in blue, green, yellow and red



Spot holder available in 85, 200 and 300 mm lengths

Base plate for attaching spot holder on which the LED spot lights can be mounted



Adapter for attaching 2-arm spots to columns and various stands (see Accessories on Page 13)





LED BACKLIGHTS



LED UNIVERSAL BACKLIGHT

This backlight (40 LEDs) can be attached to the working surface of a stereomicroscope and is characterised by daylight quality illumination, a particularly uniform illumination, high LED service life and low construction height.

The brightness of the backlight can be infinitely adjusted with the control unit.



LED NIKON BACKLIGHT

This Nikon LED backlight (equivalent to 40 LEDs) can be inserted into the bases of the Nikon SMZ 600, SMZ 800, SMZ 1000 (C-PS, C-DS stand) and SMZ 1500 (C-PS160, C-DSS, C-DSD and C-BD stands, with adapter ring) microscopes. It provides particularly bright and homogenous illumination (maximum light intensity 26,500 cd/m²) at a colour temperature of 5,500 K.

ACCESSORIES

Glass plate for flat surfaces for Petri dishes



Adapter ring for the use of the backlight for the SMZ-1500



(not illustrated:) Self-adhesive **stand feet** for raising the microscope stand

CONTROL UNITS



STANDARD CONTROL UNIT

This compact control unit has an on/off switch, an infinitely variable brightness control with "turbo range", two connection sockets for the LED lighting units.



SEGMENT CONTROL UNIT

The segment control unit also has a switch for segment control; segment selection is indicated by 4 LEDs.

- ▶ Separate button for right or left rotation of the segment
- ▶ Autorotation
- ▶ Rotation steps selectable in 45° or 90° steps (segment sizes: semi, quarter or one-eighth circles)
- ▶ Infinitely adjustable brightness
- ▶ No flicker with cameras due to direct current excitation
- ▶ Operation of two light heads possible
- ▶ Variable voltage power supply 100–240 V, ESD socket



SEGMENT CONTROL UNIT 80

This control unit enables segment selection with an on/off switch, buttons can be used to rotate the selected segments with two buttons available for direction changes. The rotary potentiometer can be used to infinitely adjust the brightness level. Compatibility is only possible with the LED ringlight 66/80. The following segments can be selected:

- ▶ Semicircle
- ▶ Quadrant
- ▶ Two opposing quadrants
- ▶ Full circle

Segment control: Different segments can be actuated by pushbutton.





HI-POWER CONTROL UNIT

Up to two Hi-Power spots can be connected to the Hi-Power control unit. It is possible to activate one spot, both spots alternatively or both spots simultaneously (the selection is shown on the display). This control unit offers infinitely variable brightness control.



COMBI CONTROL UNIT

Up to two Hi-Power spots and an LED component (e.g. a ringlight) with maximum 40 LEDs can be connected to this control unit. The LED Combi control unit can be operated in normal mode or storage mode. Any position can be entered in normal mode, while the selected positions are stored in storage mode and can be used at a later time.

ALL DEVICES – DATA AND EQUIPMENT

COMPONENTS

ILLUMINATION	LED ringlight 66/40 (619-20-002)	LED ringlight 66/80 (619-20-039)	LED ringlight 38/20 (619-20-006)
Diodes	40 white LEDs	80 white LEDs	20 white LEDs
Colour temperature	5,000 K	5,000 K	5,000 K
Dimensions	External: Ø 94 mm x 25 mm Internal: Ø 66 mm	External: Ø 112 mm x 23,5 mm Internal: Ø 66 mm	External: Ø 60 mm x 25 mm Internal: Ø 38 mm
Attachment	Circumference clamping with spring and screw		
Working distance	55–120 mm, measured from ringlight	35–120 mm or 120–300 mm, measured from ringlight	50–100 mm, measured from ringlight

ILLUMINATION	LED Hi-Power spot (619-20-037)	LED Hi-Power spot 2-arm (619-20-047)
Diodes	1 x Hi-Power LED	2 x Hi-Power LED
Dimensions	Ø 25 mm x 51 mm	
Connections	M6 thread for holder	
Attachment		Via incident light adapter

ILLUMINATION	LED spot 19 (619-20-011)	LED spot 19 2-arm (619-20-012)
Diodes	19 white LEDs	19 white LEDs per arm
Colour temperature	5,000 K	5,000 K
Dimensions	Ø 34 mm x 45 mm	Ø 34 mm x 45 mm, arm length approx. 480 mm
Attachment	With spot holder	With incident light adapter

ILLUMINATION	LED universal backlight (619-20-052)	LED Nikon backlight (619-20-036)
Diodes	40 white LEDs	LED cluster (equiv. to 40 LEDs)
Colour temperature	5,000 K	5,000 K
Dimensions	Ø 120 mm x 19 mm	Ø internal 54 mm, external 59,5 mm; external + adapter 178 mm
Slides	Frosted glass, Ø 90 mm	Flashed glass
Illuminated area	Ø 55 mm	Ø 55 mm

CONTROL UNIT	Standard (619-30-001.99)	Segment (40 LEDs) (619-30-020.99)	Segment (80 LEDs) (619-30-022.99)
Power supply	250 mA	24V DC $\pm$ 5 %, hollow plug 5.5 x 2.1 mm	500 mA
Dimensions (D x W x H)		95.5 x 64 x 29 mm	102 x 106 x 39 mm
Sum of max. individual LEDs	40	40	80
Connections		2 x 9-pin Mini-DIN-socket	D-Sub (25-pin)
On/Off switch		●	
Segment rotation pushbutton	–	2 (incl. auto-rotation)	2
Potentiometer wheel		Approx. 1–25 mA, with index holes	
ESD connection		4 mm banana plug	
Ambient temperature		10–40 °C	
Rel. humidity		30–70 %	
Mark of conformity		CE	

CONTROL UNIT	Hi-Power spot (619-30-024.99) for max. 2 Hi-Power spots	Combi (619-30-025.99) for max. 2 Hi-Power spots and 1 ringlight/backlight
Power supply	5 V DC $\pm$ 5 %, hollow plug 5.5 x 2.5 mm, 1.2 A	24 V DC $\pm$ 5 %, hollow plug 5.5 x 2.1 mm, 0,75 A
Dimensions (D x W x H)	95.5 x 64 x 29 mm	126 x 52 x 45 mm
Connections	2 x 3.5 mm jack sockets (stereo)	2 x 3.5 mm jack sockets (stereo) and 1x mini-DIN
Sum of max. individual LEDs		40
ESD connection	4 mm banana plug	4 mm banana plug
Operation	Potentiometer, pushbutton für Port 1 / Port 2 / On/Off with LED display	Potentiometer, 3 pushbuttons with LED display

ACCESSORIES

RINGLIGHT

Adapter:

66/58 for RL 66 mm on lens Ø 58 mm	595-32-258
66/60 for RL 66 mm on lens Ø 60 mm	598-32-043
66/62 for RL 66 mm on lens Ø 62 mm	598-32-127
Universal adapter for RL 66 mm on lenses Ø 38-58 mm	598-30-022
66 mm on thread M48 x 0.75 (Olympus SZ 51 and SZ 61, Nikon SMZ 445/SMZ-1)	598-32-049
66 mm on thread M55 x 0.75 (Nikon SMZ 645/660/745)	598-32-050
66/60 mm for Meiji	598-32-183
66/54 mm for Olympus SZ-III + SZ-60, SZ-1145	598-30-002
66/47 mm for Nikon (except SMZ-1)	598-30-004
66/50 mm for Olympus SD, SF, SZ with auxiliary lens	598-30-017

Diffusers:

For ringlight 66/40	619-20-040
For ringlight 66/80	619-20-044

Polarizing filter sets:

Analyser and polariser 66/40	595-20-166
Analyser and polariser 66/80	619-20-022

NIKON BACKLIGHT

Self-adhesive stand foot	619-30-050
Carrier plate	619-32-121
Glass table for stand	619-30-051

HI-POWER SPOT / HI-POWER SPOT 2-ARM

Diffuser	OPP-00006
Polariser	OPP-00007
Spot holder:	
M6/M5, L = 85 mm	619-30-041
M6/M5, L = 200 mm	619-30-042
M6/M5, L = 300 mm	619-30-043
Base plate with M5 thread	619-30-007
Adapter (for 2 arms):	
Round column adapter Ø 25, 30, 32, 35 mm	619-20-049
Olympus SZ2-ST, Nikon C-PS, C-DS	619-20-048
Leica-S-stand (for 2 arms only)	619-20-050
Nikon SMZ-1500	598-30-024

SPOT 19 / SPOT 19 2-ARM

Filter:

Blue	619-30-010
Green	619-30-011
Yellow	619-30-012
Red	619-30-013




Spot holder:

L = 85 mm	619-30-003
L = 200 mm	619-30-004
L = 300 mm	619-30-005

Base plate with M5 thread

Adapter (for 2 arms):	
Round column adapter Ø 25, 30, 32, 35 mm	619-20-049
Olympus SZ2-ST, Nikon C-PS, C-DS	619-20-048
Leica-S-stand (for 2-arms only)	619-20-050
Nikon SMZ-1500	598-30-024


BRIGHTNESS DISTRIBUTION




THE COMPREHENSIVE SYSTEM AT A GLANCE

ILLUMINATION


LED ringlight 38/20
619-20-006




LED ringlight 66/40
619-20-002



LED ringlight 66/80
619-20-039 (35–120 mm)
619-20-046 (120–300 mm)




LED Hi-Power spot
619-20-037




CONTROL UNITS


Standard (619-30-001.99)



Segment (619-30-020.99)




Combi (619-30-025.99)




Segment 80 (619-30-022.99)



Hi-Power (619-30-024.99)



Combi (619-30-025.99)




ACCESSORIES

Diffusor 66/40 (619-20-040)



Polarizing filter set 66/40
(595-20-166)




Adapter (see Page 13)



Diffusor 66/80 (619-20-044)



Polarizing filter set 66/80
(619-20-022)



Adapter (see Page 13)



Spot holder M5/M6:
85 mm (619-30-041)
200 mm (619-30-042)
300 mm (619-30-043)

Base plate M5 (619-30-007)



Diffusor (OPP-00006)
(not illustrated)
Polarizing filter (OPP-00007)
(not illustrated)

SETS

10009 **Set Ringlight 38/20**
Working dist. 50–100 mm
incl. control unit standard

10001 **Set Ringlight 66/40**
Working dist. 55–120 mm
incl. control unit standard

10049 **Set Ringlight 66/40
Segment**
Working dist. 55–120 mm
incl. control unit segment

10032 **Set Ringlight 66/80**
Working dist. 35–120 mm

10041 **Set Ringlight 66/80**
Working dist. 120–300 mm
all incl. control unit
segment 80

10035 **Set 1 Hi-Power Spot**






10037 **Set 2 Hi-Power Spots**

10043 **Set 1 Hi-Power Spot**
incl. spot holder and base plate


10044 **Set 2 Hi-Power Spots**
inkl. 2 spot holders and
2 base plates

all incl. control unit Hi-Power


(incl. specified control unit
and internal plug-in power
supply unit 100-240 V)

**LED-Hi-Power spot
2-arm 619-20-047**

**LED spot 19
619-20-011**

**LED spot 19 2-arm
619-20-012**

**LED universal backlight
619-20-052**

**LED Nikon backlight
619-20-036**



Standard (619-30-001.99)




Segment (619-30-020.99)




Combi (619-30-025.99)




Standard (619-30-001.99)



Segment (619-30-020.99)





Combi (619-30-025.99)



Adapter ring (619-32-121)




Glass plate (619-30-051)



Spot holder M5:
85 mm (619-30-003)
200 mm (619-30-004)
300 mm (619-30-005)
Base plate M5 (619-30-007)



Colour filter (see Page 13):
blue
green
yellow
red



10045 **Set Hi-Power
Spot 2-arm**
incl. control unit Hi-Power

10002 **Set 1 Spot 19**
10003 **Set 2 Spots 19**
all incl. control unit standard

10025 **Set Spot 19 2-arm**
incl. control unit standard
10027 **Set Segment
Spot 19 2-arm**
incl. control unit stegment

10056 **Set Universal
Backlight**
incl. control unit standard

10023 **Set Backlight
Nikon**
incl. control unit standard



PHOTONIC Optische Geräte GmbH & Co KG
Austria, A-1160 Vienna, Seeböckgasse 59

☎ +43-1- 486 56 91-0

📠 +43-1- 486 56 91-47

✉ sales.fo@photonic.at

www.photonic.at