
P2.6.2

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When the hot-air engine is operated as a heat engine, each engine cycle withdraws the amount of heat Q_1 from reservoir 1, generates the mechanical work W and transfers the difference $Q_2 = Q_1 - W$ to reservoir 2. The hot-air engine can also be made to function as a refrigerating machine while operated in the same rotational direction by externally applying the mechanical work W . In both cases, the work W_F converted into heat in each cycle through the friction of the piston in the cylinder must be taken into consideration.

In order to determine the work of friction W_F in the experiment P2.6.2.1, the temperature increase DT_F in the cooling water is measured while the hot-air engine is driven using an electric motor and the cylinder head is open.

The experiment P2.6.2.2 determines the efficiency

$$\eta = \frac{W}{W + Q_2}$$

of the hot-air engine as a heat engine. The mechanical work W exerted on the axle in each cycle can be calculated using the external torque N of a dynamometrical brake which brakes the hot-air engine to a speed f . The amount of heat Q_2 given off corresponds to a temperature increase DT in the cooling water.

The experiment P2.6.2.3 determines the efficiency

$$\eta = \frac{Q_2}{Q_1 - Q_2}$$

of the hot-air engine as a refrigerating machine. Here, the hot-air engine with closed cylinder head is driven using an electric motor and Q_1 is determined as the electrical heating energy required to maintain the cylinder head at the ambient temperature

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
*dodatkowo zalecane

Thermodynamic cycles are often described as a closed curve in a pV diagram (p: pressure, V: volume). The work added to or withdrawn from the system (depending on the direction of rotation) corresponds to the area enclosed by the curve.

In the experiment P2.6.2.4, the pV diagram of the hot air engine as a heat engine is recorded using the computer-assisted measured value recording system CASSY. The pressure sensor measures the pressure p in the cylinder and a displacement sensor measures the position s, from which the volume is calculated, as a function of the time t. The measured values are displayed on the screen directly in a pV diagram. In the further evaluation, the mechanical work performed as piston friction per cycle

$$W = -\int p \cdot dV$$

and from this the mechanical power

$$P = W \cdot f$$

f : no-load speed

are calculated and plotted in a graph as a function of the no-load speed.