


SAUTER TI-AC

Lever operated test stand for hardness testing with base plate made of glass

KERN

Product information

GTIN/EAN number	4045761140233
Trade number	90248000
REACH classification	no prohibited substances

Pictograms

STANDARD



Date of creation: 20.06.2026

Category

Brand	Sauter
Product category	Hardness testing
Product group	Manual Shore test stand
Product family	TI

Measuring System

Hardness scales	Shore A
Test force Hardness measurement	10 N

Construction

Dimensions completely mounted (W×D×H)	150×200×330 mm
Base plate diameter (Ø)	75 mm
Material housing	stainless steel
Test object height [Max]	60 mm
Traverse distance with lever movement / release [Max]	15 mm

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	380×255×195 mm
Shipping method	Parcel service
Net weight approx.	7 kg
Gross weight approx.	8 kg
Shipping weight	7,5 kg
Shipping volume weight	3,8 kg