


SAUTER TI-D

Lever operated test stand for hardness testing with base plate made of glass

KERN

Product information

GTIN/EAN number	4045761140240
Trade number	90248000
REACH classification	no prohibited substances

Pictograms

STANDARD



Date of creation: 20.06.2026

Category

Brand	Sauter
Product category	Hardness testing
Product group	Manual Shore test stand
Product family	TI

Measuring System

Hardness scales	Shore D
Test force Hardness measurement	50 N

Construction

Dimensions completely mounted (W×D×H)	150×200×400 mm
Base plate diameter (Ø)	75 mm
Material housing	stainless steel
Test object height [Max]	60 mm
Traverse distance with lever movement / release [Max]	15 mm

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	355×255×180 mm
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	10 kg
Shipping weight	9,5 kg
Shipping volume weight	3,3 kg