

KERN OZB-A4571



Professional illumination units guarantee outstanding, uniform and strong illumination



Date of creation: 11.07.2026

Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IR

Construction

Ring illumination diameter	6,1 cm
Beleuchtungsart	LED
Illuminance	4W
Colour temperature	7000-11000 K
Illumination dimmable	✓

Power Supply

Supplied power supply	Power supply unit
	EURO
	AUS
Plug-in power supply / adapter for countries - optional	UK
	US
	CH

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Approval

CE mark	✓
---------	---

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	235×200×50 mm
Shipping method	Parcel service
Net weight approx.	0,45 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg
Shipping volume weight	0,5 kg

Product information

ATTR_warrantyTime	24 mon
Trade number	94054239
GTIN/EAN number	4045761218420

Pictograms

STANDARD

