

SAUTER TI-HED

Test stand for hardness testing Shore A and D

KERN

Trade number	90248000
REACH classification	no prohibited substances

Pictograms

STANDARD



Date of creation: 22.06.2026

Category

Brand	Sauter
Product category	Hardness testing
Product group	Manual Shore test stand
Product family	TI-HE

Measuring System

Hardness scales	Shore D
Test force Hardness measurement	50 N

Construction

Dimension (W×D×H)	200×200×470 mm
Dimension housing (W×D×H)	200×200×470 mm
Dimensions completely mounted (W×D×H)	200×200×470 mm
Base plate diameter (Ø)	115 mm
Material housing	steel
Type of drive	manual crank handle
Test object height [Max]	50 mm
Traverse distance with lever movement / release [Max]	20 mm

Packing & Shipping

Delivery time	1 d
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	13 kg
Shipping weight	12,35 kg

Product information

GTIN/EAN number	4045761274044
-----------------	---------------